



DRAGON MEDICAL ONE / DRAGON PRO ANYWHERE HOW TO ADD NEW WORDS AND PHRASES

1. Click the **DragonBar menu > Manage Vocabulary...**:



2. The Manage Vocabulary screen will be displayed; click the **PLUS SIGN (+)** to add a new word:



3. Type the word exactly as you would like it to appear when recognised (with appropriate capitalisation, etc), for example "VoiceX"

4. Click one of the following buttons to add the word:

Dragon Medical One

Back Manage Vocabulary

Add Word

Enter the word you want to add:

VoiceX

a Use Default Pronunciation

Dragon provides the default pronunciation based on the word's spelling. If you pronounce it differently from the spelling you can click the microphone icon here and say the word as you would normally pronounce it.

b

c Enter the word as you would pronounce it

- a. Use Default Pronunciation:** Dragon will estimate the word's pronunciation – this is sufficient for simple words that are phonetic and do not include any special characters (eg. "Richmond")
- b. Microphone Icon:** Dragon will ask you to say the word to train it – this is a better option if the word isn't strictly phonetic, but still sounds like it looks for the most part (eg. "Caulfield")
- c. Enter the word as you would pronounce it:** Enter a spoken form to clearly define how you pronounce the word – this is the best option for words that contain special characters or are simply not dictated as they are spelled. Use existing words where possible and upper case characters to indicate dictation of a letter:

Dragon Medical One

Back Manage Vocabulary

Add Word

Enter the word you want to add:

VoiceX

Pronunciation:


voice X

If how a word is pronounced is different from how it is written, you should also enter a pronunciation. For example:


Word: "IEEE"
Pronunciation: "I triple E"

✓ Save

5. Once a word has been added, you can provide additional training (if required) by returning to **Manage Vocabulary** and selecting the word, following by **Train**:



6. You can create multiple-word phrases by following the same procedure as individual words. Lengthy phrases, or those that are tricky to dictate (due to symbols, etc), can be shortened using option **c** in **Step 4**:



7. Custom words can be deleted by returning to **Manage Vocabulary** and selecting the word(s), following by the **Trash Can** icon:

